

ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

Factory address: No 1, 2 Factory, No. 5, Hutang Road, Xiaya Town, Jiande City, Hangzhou, Zhejiang Province, China

Office address: No 1, 2 Factory, No. 5, Hutang Road, Xiaya Town, Jiande City, Hangzhou, Zhejiang Province, China

TEL: 0086-571-86676309

E-mail: info@geyi-medical.com

www.geyi-medical.com



STAY TUNED

Product upgrade coming soon



Single-use

Digital Flexible Ureteroscope




Single-use
Digital Flexible Ureteroscope

3.6Fr instrument channel, compatible with standard urology instruments.



- 1:1 coaxial steering, more smooth to operate.
- Two shortcut keys, support user customization.

type-c plug,
electrical connection is more stable and reliable.



275°Bidirectional bending




Probe with smooth bevel design
smooth insertion and operation



Product
Feature

1. 275° bidirectional bending, 1:1 coaxial steering.
2. Handle keys can be preset with up to 15 functions.
3. Probe with smooth bevel design, smooth insertion and operation.
4. Front LED light source, High sensitivity CMOS, excellent color fidelity.
5. 160,000 pixel imaging

Hd camera restores ultimate color



Patented design with 15° bevel Angle

Compatible with multiple hosts



Product
Parameter

Model	Max-insert Diameter (mm)	Bending angle	Working length (mm)	Working channel (mm)	Viewing angle	Depth of field (mm)
GY-UR9.3	9.3Fr(3.1mm)	275°	670mm	3.6Fr(1.2mm)	120°	3-50mm
GY-UR8.4	8.4Fr(2.8mm)	275°	670mm	3.6Fr(1.2mm)	120°	3-50mm
GY-UR7.5	7.5Fr(2.5mm)	275°	670mm	3.4Fr(1.1mm)	120°	3-50mm
GY-UR6.3	6.3Fr(2.1mm)	275°	670mm	3.4Fr(1.1mm)	120°	3-50mm